

Scar|Fx®

Product Description: Scar Fx® Silicone Sheeting is made of 100% healthcare grade silicone, manufactured via a proprietary method allowing the sheeting to be flexible, durable, and have a natural tack without the addition of adhesives.

Product Indications: Scar Fx® Silicone Sheeting is used to improve the appearance of surgical and traumatic scars (including normal post-op scars, hypertrophic scars, and keloid scars) and burns. Scar Fx® is applied to a clean and dry scar, for up to 8-12 hours nightly, and can be secured to the skin with medical tape as needed. Scar Fx® creates an occlusive barrier promoting skin moisture. This barrier supports the mechanism of action, which includes hydration and occlusion.

SECC Chart

SAFETY	<p>Scar Fx® Silicone Sheeting is classified by the FDA/CE as an OTC Class 1 Medical Device, which is 510K/PMA-Exempt.</p> <p>Scar Fx® is indicated for use on epithelialized (closed) skin. Scar Fx® has been in the market since 1998, is sold worldwide, and has an excellent record of safety with no adverse events reported. Silicone is naturally hypoallergenic, making this product safe for use on all skin types. Scar Fx® is not recommended for patients under the age of three due to the risk of choking hazard. Please consult your physician.</p>
EFFICACY	<p>The International Advisory Panel on Scar Management recommends silicone sheeting as the "... standard care for plastic surgeons."⁽¹⁾ for non-invasive scar management.</p> <p>Silicone Sheeting has been studied in excess of thirty years, and is considered the first line of defense against abnormal scarring. Studies conclude that "It is much more efficient to prevent hypertrophic scars than to treat them."⁽¹⁾</p> <p>Scar Fx® Silicone Sheeting is made with the highest quality ingredients and to the highest specifications to achieve the maximum result for patients seeking to manage new and old scars.</p>
COST EFFECTIVE	<p>Scar Fx® provides a tremendous value to both the physician and patient based on its quality, durability, and ease of use. A single Scar Fx® silicone sheet can be used for the duration of treatment (30-45 days for most normal surgical scars).</p> <p>Physicians receive preferential pricing due to the positioning of the product as physician-recommended. Scar Fx® Silicone Sheeting has an industry-low return rate under 0.05% compared to a 2% industry average.</p>
COMPLIANCE	<p>Scar Fx® Silicone Sheeting's regimen was designed to maximize compliance. Scar Fx® is applied in the evening and worn while sleeping. In the morning, the sheeting is removed and the process is repeated.</p>

(1) Mustoe TA, Cooter RD, Gold MH, Hobbs FD, Ramelet AA, Shakespeare PG, Stella M, Téot L, Wood FM, Ziegler UE; International Advisory Panel on Scar Management. Plast Reconstr Surg. 2002 Aug;110(2):560-71.

Manufactured in the USA for:

Scar Heal
13191 Starkey Road, Bldg. 11
Largo, Florida 33773-1438



Document	Revision
# SH-SSS-SECC	03

